

MILLENSYS



## **3D Secondary Capture SOP Class**

© Copyright by MILLENSYS

## 1. Introduction

The Vtools Family of applications will create functional images from special applications. Those will be encoded as Standard extended SOP Class. Please see the following tables for a selected overview of supplied Standard and Private attributes.

SC SOP Class		
Attribute Name	Tag	Value
Specific Character Set	(0008,0005)	<as defined with Study>
Image Type	(0008,0008)	'ORIGINAL\SECONDARY\OTHER'
Instance Creation Date	(0008,0012)	
Instance Creation Time	(0008,0013)	
Instance Creator UID	(0008,0014)	
SOP Class UID	(0008,0016)	<depending on the modality type>
SOP Instance UID	(0008,0018)	<new UID>
Study Date	(0008,0020)	<yyyymmdd><as defined with Study>
Series Date	(0008,0021)	<yyyymmdd><as defined with Study>
Acquisition Date	(0008,0022)	<yyyymmdd><file creation date>
Content Date	(0008,0023)	<yyyymmdd> <file creation date>
Acquisition Date time	(0008,002A)	<yyyymmddhhmmss><file creation date and time>
Study Time	(0008,0030)	<hhmmss ><as defined with Study>
Series Time	(0008,0031)	<hhmmss ><as defined with Study>
Acquisition Time	(0008,0032)	<hhmmss ><file creation date>
Content Time	(0008,0033)	<hhmmss ><file creation date>
Accession Number	(0008,0050)	<as defined with Study>
Modality	(0008,0060)	CT or MR or US
Conversion Type	(0008,0064)	DV
Manufacturer	(0008,0070)	'MillenSys Systems'
Institution Name	(0008,0080)	<Defined in the site configuration>
Institution Address	(0008,0081)	<Defined in the site configuration>
Referring Physician's Name	(0008,0090)	<as defined with Study>
Station Name	(0008,1010)	
Study Description	(0008,1030)	<as defined with Study>
Series Description	(0008,103E)	'Series Created VT3dView, Volume Captured Images'

<b>Performing Physicians' Name</b>	(0008,1050)	
<b>Admitting Diagnoses Description</b>	(0008,1080)	
<b>Manufacturer's Model Name</b>	(0008,1090)	'Vision Tools '+software verison
<b>Patient's Name</b>	(0010,0010)	<as defined with Patient>
<b>Patient ID</b>	(0010,0020)	<as defined with Patient>
<b>Patient's Birth Date</b>	(0010,0030)	<as defined with Patient>
<b>Patient's Birth Time</b>	(0010,0032)	
<b>Patient's Sex</b>	(0010,0040)	<as defined with Patient>
<b>Other Patient IDs</b>	(0010,1000)	
<b>Other Patient Names</b>	(0010,1001)	
<b>Patient's Age</b>	(0010,1010)	
<b>Patient's Size</b>	(0010,1020)	<as defined with Study>
<b>Ethnic Group</b>	(0010,2160)	
<b>Occupation</b>	(0010,2180)	
<b>Additional Patient's History</b>	(0010,21B0)	
<b>Patient Comments</b>	(0010,4000)	<as defined with Patient>
<b>Body Part Examined</b>	(0018,0015)	<as defined with Series>
<b>Study Instance UID</b>	(0020,000D)	<as defined with Study>
<b>Series Instance UID</b>	(0020,000E)	<new UID>
<b>Study ID</b>	(0020,0010)	<as defined with Study>
<b>Series Number</b>	(0020,0011)	<highest series number + 1>
<b>Instance Number</b>	(0020,0013)	<highest image number + 1>
<b>Image Comments</b>	(0020,4000)	'Processed Image with Vision Tools ' +Software version
<b>Samples per Pixel</b>	(0028,0002)	1
<b>Photometric Interpretation</b>	(0028,0004)	MONOCHROME2 or RGB
<b>Rows</b>	(0028,0010)	
<b>Columns</b>	(0028,0011)	
<b>Pixel Aspect Ratio</b>	(0028,0034)	-
<b>Bits Allocated</b>	(0028,0100)	16 bit or 24 bit
<b>Bits Stored</b>	(0028,0101)	8-12 bit
<b>High Bit</b>	(0028,0102)	7-11
<b>Pixel Representation</b>	(0028,0103)	0
<b>Smallest Image Pixel Value</b>	(0028,0106)	
<b>Largest Image Pixel Value</b>	(0028,0107)	
<b>Window Center</b>	(0028,1050)	
<b>Window Width</b>	(0028,1051)	



<b>Rescale Intercept</b>	(0028,1052)	
<b>Rescale Slope</b>	(0028,1053)	
<b>Lossy Image Compression</b>	(0028,2110)	
<b>Pixel Data</b>	(7FE0,0010)	binary

MILLENSYS