


# *Vision Tools*

## *CTView*

**Complete Digital Station for CT Machines**



**Huge Archive and Rapid Operations for MultiSlice CTs**  
**CD/DVD publishing with dicom viewer**  
**Visual Film printing on DICOM and A4 printers**  
**Orthogonal, Oblique, Curved MPR/MIP Option**  
**3D powerful reconstruction tools option**  
**Visual Hanging protocols and Multi-Monitor Support**  
**Connect non-dicom CT to PACS to send and print dicom**  
**Support Video, Digital and Dicom CTs**

**Display Scan locators and Pilot images**  
**Acquisition and printing with remote control**  
**Image processing, annotations and CT measurements**  
**Reporting with voice dictation**  
**DICOM Send,QR,Print, WL/MPPS Options**  
**Send/Receive Studies via E-mail Option**

# Vision Tools CTView

## Connectivity

Accepts video and DICOM studies from single or multiple CT machines. Supports ACR/NEMA studies using MiXLink and MiTOS gateways.

## Images

Huge number of image processing, filters, annotations, measurements, window adjustment and display functions. MPR / 3D reconstruction option.

## Reporting

Fully functional report editor with template and structured reporting support. voice notes recording during and after exam.

## DICOM

Built in autorun DICOM viewer featuring all needed functions for viewing. Supports CT modality and SC Export DICOM Archive, DICOM send/receive, DICOM Q/R, Basic Modality worklist, MPPS and DICOM print, ... etc.

## Patients


Detailed patient database with quick and advanced search functions. Archiving of studies on any archiving media.

## Printing

Supports multi format and film queuing. Supports dry/wet laser imagers and windows compatible printers. visual filming for quickest film printing with online processing over the entire film.

## Options

Remote control for video Interface. Multi Monitor support. Multilingual Support. Compatible with PACS environment. PACS/RIS integration with MiRIS. Vision Tools 3DView for 3D reconstruction. MiXLink or MiTOS for ACR/NEMA support.



millensys reserves the rights to modify the design and specifications contained herein without prior notice. Please refer to millensys web site for the most current information